

“Jai Sri Gurudev”
Sri Adhichunchanagiri Shikshana Trust ®


B.G.S.INSTITUTE OF TECHNOLOGY

B.G NAGARA - 5 71448
MANDYA



**ENGINEERING CHEMSITRY
LABORATORY MANUAL**

I/II SEMESTER B.E.



**PREPARED BY :
FACULTY, DEPARTMENT OF CHEMISTRY**

Engineering Chemistry Laboratory

Sub Code: 18CHEL16/18CHEL 26

Time : 03 Hours

Marks : 80

(For Examination, any one experiment from Part- A and a common experiment from Part - B shall be set)

PART – A

| SL. No. | Experiment | Page No. |
|---------|--|----------|
| 01 | Potentiometric estimation of FAS using standard $K_2Cr_2O_7$ solution | 4 |
| 02 | Colorimetric estimation of Copper | 8 |
| 03 | Conductometric estimation of an Acid mixture using standard NaOH solution | 12 |
| 04 | Determination of pKa of a weak acid using pH Meter | 16 |
| 05 | Determination of Viscosity co-efficient of a given liquid | 20 |
| 06 | Flame Photometric estimation of Sodium and Potassium in the given sample of Water. | 24 |

PART – B

| SL. No. | Experiment | Page No. |
|---------|---|----------|
| 01 | Determination of Total Hardness of a sample of Water using Disodium salt of EDTA. | 26 |
| 02 | Determination of CaO in the given sample of Cement solution by Rapid EDTA method. | 30 |
| 03 | Determination of Percentage of Copper in Brass using standard Sodium thiosulphate solution | 34 |
| 04 | Determination of Iron in the given sample of Haematite ore solution using $K_2Cr_2O_7$ Crystals by external indicator method. | 38 |
| 05 | Determination of COD of the given Industrial Waste water sample. | 42 |
| 06 | Determination of Total Alkalinity of a given Water Sample using standard Hydrochloric acid. | 48 |

References Books:

1. Laboratory manual in Engineering Chemistry Sudharani, Dhanpatrai Publishing Company.
2. Vogel's Text Book of Quantitative Chemical Analysis revised by G.H.Jeffery, J. Bassett, J. Mendham and R.C Denney.

3. INSTRUCTIONS FOR INSTRUMENTAL METHOD OF ANALYSIS

- ❖ Read the procedure before commencement of experiment and do the experiment according to the procedure.
- ❖ Instruments used for Part – A experiments are:
- ❖ pH meter, Potentiometer, Colorimeter, Conductivity meter and Viscometers.
- ❖ Wash the micro burette, pipette and beaker with tap water and rinse with distilled water before commencement of the experiments.
- ❖ Wash the electrodes or cells with distilled water and wipe them with tissue paper before commencement of the experiments.
- ❖ Switch off the instruments after experiment is completed.
- ❖ Wash the electrodes or cells with distilled water and dip them in distilled water taken in a 100 cm³ beaker after the experiments is completed.
- ❖ Wash the micro burette, pipette, standard flask and beaker with tap water again after the experiment is completed and return them to the instructor.

INSTRUCTIONS FOR VOLUMETRIC METHOD OF ANALYSIS

- ❖ Read the procedure before commencement of the experiment and do the experiment according to the procedure.
- ❖ Handle the glass wares with care.
- ❖ Apparatus used: burette, pipette, standard flask, conical flask, beaker, funnel, washing bottle.
- ❖ Wash the all apparatus with tap water and rinse them with distilled water before commencement of the experiment.
- ❖ Wash all the apparatus with tap water after the experiment is completed and return them to the instructor

PART-A

INSTRUMENTAL
METHOD OF ANALYSIS

| Volume of $\text{K}_2\text{Cr}_2\text{O}_7$ | E.M.F (E) in mV | ΔE | ΔV | $\Delta E/\Delta V$ |
|--|----------------------------|------------|------------|---------------------|
| 0.0 | | ** | 0.5 | |
| 0.5 | | | | |
| 1.0 | | | | |
| 1.5 | | | | |
| 2.0 | | | | |
| 2.5 | | | | |
| 3.0 | | | | |
| 3.5 | | | | |
| 4.0 | | | | |
| 4.5 | | | | |
| 5.0 | | | | |
| 5.5 | | | | |
| 6.0 | | | | |
| 6.5 | | | | |
| 7.0 | | | | |
| 7.5 | | | | |
| 8.0 | | | | |
| 8.5 | | | | |
| 9.0 | | | | |
| 9.5 | | | | |
| 10.0 | | | | |

INSTRUMENTAL METHODS OF CHEMICAL ANALYSIS

PART-A

EXPERIMENT NO. 1

Date:

| | | | | | | | |
|--|--|--|--|--|--|--|--|
| | | | | | | | |
|--|--|--|--|--|--|--|--|

POTENTIOMETRIC ESTIMATION OF FAS USING STANDARD $K_2Cr_2O_7$ SOLUTION.

THEORY:

When the electrode is dipped in the solution, the potential of an electrode depends on the concentration of the ions in the electrolyte. Therefore determination of maximum potential can serve as an indicator in potentiometric analysis. Determination of change of potential with concentration gives more accurate results than volumetric analysis.

The given electrode is combined with saturated calomel electrode and dipped in the electrolyte, an electro chemical cell is formed. The emf of this cell is measured using potentiometer. During the titration, the change in emf of the cell or concentration of ions in the electrolyte can be measured. Then $\Delta E/\Delta V$ is plotted against the volume of the titrant, there will be a peak corresponding to the end point of the reaction.

\therefore emf of the cell is calculated using Nernst equation,


$$E_{\text{cell}} = E^{\circ}_{\text{cell}} + \frac{0.0591}{n} \log \frac{[\text{Pr oducts}]}{[\text{Re ac tan ts}]}$$

Where E_{cell} = emf of the cell, E°_{cell} standard emf of the cell,

n = Number of electrons involved in the reaction,

[Products] = concentration of the products and

[Reactants] = concentration of the reactants.



Scale :

X-axis, 1.0 cm = _____ cm^3

Y-axis, 1.0 cm = _____ units

\therefore Volume of $\text{K}_2\text{Cr}_2\text{O}_7$ required = $V_{\text{K}_2\text{Cr}_2\text{O}_7} = \text{_____ cm}^3$

Note:

1. Normality of $\text{K}_2\text{Cr}_2\text{O}_7 = N_{\text{K}_2\text{Cr}_2\text{O}_7} = 0.1 \text{ N}$

2. Equivalent weight of Mohr's salt = 392

\therefore Normality of Mohr's salt solution = $\frac{N_{\text{K}_2\text{Cr}_2\text{O}_7} \times V_{\text{K}_2\text{Cr}_2\text{O}_7}}{V_{\text{Mohr's salt}}}$

$N_{\text{Mohr's salt}} = \frac{\text{_____}}{25} = \text{_____ N}$

Weight per liter of

Mohr's salt crystals = $X = N_{\text{Mohr's Salt}} \times \text{equivalent weight of Mohr's salt}$

$X = \text{_____} \times 392 = \text{_____ g}$


\therefore Amount of Mohr's salt present in the given solution is =g.

RESULT: Amount of Mohr's salt present in the given solution isg.

PROCEDURE:

Transfer the given amount of Mohr's salt solution (25 cm³) in to a clean 100 cm³ beaker. to that add one test tube of dilute sulphuric acid.

Switch on the potentiometer and adjust the temperature knob to the lab temperature. Connect the platinum and calomel electrodes to the potentiometer. Wash the electrodes with distilled water and immerse the electrodes in to the Mohr's salt solution taken in a beaker. Measure the emf of the solution in terms of **mV**. Add 0.5 cm³ of potassium dichromate solution at a time taken in a micro burette into the Mohr's salt solution. Stir the solution and measure the emf of the solution and continue the addition potassium dichromate added until the sudden raise in the emf. Then take few more readings and take the difference in the emf of solution (ΔE) and the volume of potassium dichromate added (ΔV). Plot the graph, $\Delta E/\Delta V$ against the volume of potassium dichromate added. From the graph, determine the volume of potassium dichromate required to rapid increase in the emf of the solution. Calculate the amount of Mohr's salt given in the solution.



(FAS)

RESULT: Amount of Mohr's salt present in the given solution is
.....g.

OBSERVATIONS AND CALCULATIONS

| Volume of copper sulphate in cm^3 | Concentration of copper sulphate in mg / cm^3 | Concentration of copper in mg / cm^3 | Absorbance (A) or OD at 620nm |
|--|---|--|-------------------------------|
| 5 | $3.93 \times 5 = 19.65$ | $05 \times 1 = 05$ | |
| 10 | $3.93 \times 10 = 39.30$ | $10 \times 1 = 10$ | |
| 15 | $3.93 \times 15 = 58.95$ | $15 \times 1 = 15$ | |
| 20 | $3.93 \times 20 = 78.60$ | $20 \times 1 = 20$ | |
| 25 | $3.93 \times 25 = 98.25$ | $25 \times 1 = 25$ | |
| Test solution | | | |


Note: Molecular weight of copper sulphate pentahydrate = 249.54.

1. Atomic weight of copper = 63.54

1.0 cm^3 of stock solution = 3.93 mg of copper sulphate

249.54 g of $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ = 63.54 g of copper

$\therefore 3.93 \text{ mg of } \text{CuSO}_4 = \frac{63.54 \times 3.93}{249.54} = 1.0 \text{ mg of Copper.}$



Scale :

X-axis, 1.0 cm = _____ mg/cm^3

Y-axis, 1.0 cm = OD

EXPERIMENT NO. 2Date:

| | | | | | | | |
|--|--|--|--|--|--|--|--|
| | | | | | | | |
|--|--|--|--|--|--|--|--|

DETERMINATION OF COPPER BY COLORIMETRIC METHOD**THEORY:**

When a monochromatic light of intensity I_o is incident on a transparent medium, a part of light is absorbed by the media (I_a), a part of light is reflected (I_r) and the remaining part of light is transmitted (I_t).

$$\therefore I_o = I_a + I_t + I_r \quad (\text{for a glass - air interface, } I_r \text{ is negligible})$$

$$\therefore I_o = I_a + I_t$$

The ratio of $I_t/I_o = T$ called transmittance and $\log I_o/I_t = A$ which is called absorbance or optical density. The relationship between absorbance (A), concentration (c) (expressed in mol/dm^3) and path length (t) (expressed in cm) is given by Beer - Lambert's law.

$$A = \log I_o / I_t = \epsilon ct.$$

Where ϵ is the molar extension coefficient which is a constant for a given substance at a given wave length and t is the path length of the cell tube which is also constant for given cell tube. Therefore, **$A \propto c$** . Hence plot the graph, absorbance against concentration of the solution gives a straight line.

A series of standard solution of copper salt and test solution is treated with ammonia to get blue cuprammonium complex and is diluted to a definite volume. The absorbance of each of these solution and test solution is measured at **620 nm** since the complex shows maximum absorbance at this wavelength. The absorbance values are plotted against concentration of copper get a calibration curve and concentration of copper in test solution is determined from the graph.

$\therefore 1.0 \text{ cm}^3 \text{ of CuSO}_4 = 3.93 \text{ mg of CuSO}_4 = 1.0 \text{ mg of copper}$

C = concentration of copper in test solution = _____ mg / cm^3

V cm^3 of copper sulphate = C mg of copper


$$\therefore V \text{ cm}^3 \text{ of copper sulphate} = V = \frac{C}{1.0} = \frac{\quad}{1.0} = \quad \text{cm}^3$$

Hence, the volume of test solution given = V = _____ cm^3

PROCEDURE:

Transfer the given copper sulphate solution (Stock solution) to a burette and draw out 5, 10, 15, 20 and 25 cm³ of the solution into 50 cm³ standard flasks. Add 5 cm³ of 1:1 ammonia solution to each of them and also to the test solution given in a separate 50 cm³ standard flask. Make up the solution to the mark by adding distilled water and mix well for uniform concentration. Prepare a blank solution by adding 5 cm³ of 1:1 ammonia solution into a 50 cm³ standard flask make up the solution to the mark with distilled water and mix well for uniform concentration.

Switch on the colorimeter and warm it for 5 minutes. Using the filter knob (B), select the wave length or $\lambda_{\max}=620$ nm for copper sulphate solution. Take the blank solution in the cell tube and keep it in the colorimeter and adjust zero absorbance for 100% transmittance using the knob (A). Then measure the optical density or absorbance for each standard solution and test solution using a photo electric colorimeter. A straight line graph is obtained by plotting optical density or absorbance against the concentration of copper. From the graph, determine the concentration of copper in the test solution and calculate the volume of copper sulphate in the given test solution.



Preparation of blank solution: - Blank solution can be prepared by adding 5.0 cm³ of 1:1 ammonia and distilled water in 50 cm³ standard flask. Mix well for uniform concentration.


Results:

1. The volume of test solution given = _____ cm³
2. Absorbance of test solution = _____
3. Concentration of copper in test solution = _____ mg/ cm³
4. Concentration of copper sulphate solution = _____ g/cm³

OBSERVATIONS AND CALCULATIONS

Estimation of acid mixture.

| Vol. of NaOH in cm^3 | Conductance $\text{Ohm}^{-1} \text{cm}^{-1}$ | Vol. of NaOH in cm^3 | Conductance $\text{Ohm}^{-1} \text{cm}^{-1}$ |
|----------------------------------|---|----------------------------------|---|
| 0.0 | | 5.5 | |
| 0.5 | | 6.0 | |
| 1.0 | | 6.5 | |
| 1.5 | | 7.0 | |
| 2.0 | | 7.5 | |
| 2.5 | | 8.0 | |
| 3.0 | | 8.5 | |
| 3.5 | | 9.0 | |
| 4.0 | | 9.5 | |
| 4.5 | | 10.0 | |
| 5.0 | | 10.5 | |



EXPERIMENT NO. 3Date:

| | | | | | | | |
|--|--|--|--|--|--|--|--|
| | | | | | | | |
|--|--|--|--|--|--|--|--|

**CONDUCTOMETRIC ESTIMATION OF ACID MIXTURE USING STANDARD
SODIUM HYDROXIDE****THEORY:**

Measurement of conductance can be employed to determine the end point in acid-base titrations. In conductmetric titrations, there is a sudden change in conductance of the solution near the end point. The end point is determined graphically by plotting conductance (ordinate) against sodium hydroxide (abscissa).

Conductometric titration may be applied for the determination of acid & acids present in a mixture. In the titration of a mixture of acid with a strong base (sodium hydroxide), the conductance falls until the strong acid is neutralized completely, and the weak acid undergo dissociate and gets neutralized. This results in the increases in conductance of the solution as the weak acid is consumed and converted into salt. When the neutralization of second acid is complete, there is a steep increase in conductance due to ions provided by the strong base.

PROCEDURE:

Transfer the given acid mixture (25 cm³) into a clean 100 cm³ beaker. Fill standard NaOH (0.2N) solution into 10 cm³ burette. Place the conductivity cell inside the beaker, so that the electrodes are immersed completely in the acid solution. Connect the cell to the conductometer and measure the conductance. Add 0.5 cm³ of NaOH each time, till there is a rise in conductance. After each addition stir the solution and note down the corresponding decrease in the conductivity of the solution. As the titration proceeds the conductance is gradually a decreases and then increases slowly and finally increases sharply. Take 5-6 readings when conductance begins to increase sharply. Plot a graph of conductance against volume of NaOH added.

Estimation of HCl

Volume of NaOH required to neutralize HCl = V_1 =cm³

$$\text{Normality of HCl} = \frac{N_{\text{NaOH}} \times V_1}{25}$$

$$= \frac{0.2 \times \dots\dots\dots}{25} = \dots\dots\dots \text{N}$$

Weight per liter of HCl = X = N_{HCl} X Equivalent weight of HCl

$$= \dots\dots\dots \text{X } 36.5$$

$$X = \dots\dots\dots \text{g.}$$

Estimation of CH₃COOH

Volume of NaOH required to neutralize CH₃COOH = ($V_2 - V_1$) =cm³

$$\text{Normality of CH}_3\text{COOH} = \frac{N_{\text{NaOH}} \times V_{(V_2 - V_1)}}{25}$$



$$= \frac{0.2 \times \dots\dots\dots}{25} = \dots\dots\dots \text{N}$$

Weight per liter of CH₃COOH = X = $N_{\text{CH}_3\text{COOH}}$ X Equivalent weight of CH₃COOH

$$= \dots\dots\dots \text{X } 60.05$$

$$X = \dots\dots\dots \text{g.}$$

From the graph determine the volume of NaOH required be V_1 cm³ for strong acid. As the titration proceeds, the second neutralization point V_2 will be obtained. The Difference $V_2 - V_1$ will give the amount of base neutralization for weak acid. Calculate the weight per liter of HCl and CH₃COOH.



RESULT: a) The weight of HCl/dm³ _____g
b) The weight of CH₃ COOH/dm³ _____g

OBSERVATIONS AND CALCULATIONS

| Volume of NaOH (cm ³) | pH | ΔP^H | ΔV | $\Delta P^H / \Delta V$ |
|-----------------------------------|----|--------------|------------|-------------------------|
| 0.0 | | | 0.5 | |
| 0.5 | | | | |
| 1.0 | | | | |
| 1.5 | | | | |
| 2.0 | | | | |
| 2.5 | | | | |
| 3.0 | | | | |
| 3.5 | | | | |
| 4.0 | | | | |
| 4.5 | | | | |
| 5.0 | | | | |
| 5.5 | | | | |
| 6.0 | | | | |
| 6.5 | | | | |
| 7.0 | | | | |
| 7.5 | | | | |
| 8.0 | | | | |
| 8.5 | | | | |
| 9.0 | | | | |
| 9.5 | | | | |
| 10.0 | | | | |

RESULTS:

1. Equivalence Point = _____
2. Half E.P = _____
3. pH at E.P = _____
4. pH at H.E.P. = _____
5. pKa of the given acid = _____


EXPERIMENT NO. 4Date:

| | | | | | | | |
|--|--|--|--|--|--|--|--|
| | | | | | | | |
|--|--|--|--|--|--|--|--|

DETERMINATION OF pK_a VALUE OF A WEAK ACID USING pH METER**THEORY:**

The strength of an acid is experimentally measured by determining its equilibrium constant or dissociation constant (K_a). Strong acids are completely ionized in aqueous solution and hence it is impossible to determine the dissociation constant of strong electrolytes. But in case of weak electrolytes are partially ionized in aqueous solution and hence the equilibrium constant of weak electrolytes can be experimentally determined.

Let us consider weak acid like formic acid, which ionizes, freely in aqueous solution as,





pK_a is a modern method of expressing the strength of weak acid ,

ie., pK_a = - log₁₀ K_a

pK_a value is determined by measuring the changes in pH of acid solutions at different amount of the base added. In the titration of an acid with a base, the pH of the solution rises gradually at first and then more rapidly. At the equivalence point there is a very sharp increase in pH for a very small quantity of base added. Beyond the equivalence point, the pH increases slightly on addition of base. The titration curve is obtained by plotting changes in pH against different volume of the base added and from the graph, we obtain the equivalence point of an acid. The pK_a value of an acid can be calculated using Henderson equation,

$$\text{pH} = \text{pK}_a + \log_{10} \frac{[\text{Salt}]}{[\text{Acid}]}$$

At half equivalence point [salt] = [acid] and hence pH at half equivalence point gives the pK_a value of weak acid, i.e. pH = pK_a .



PROCEDURE:

Transfer given quantity of weak acid (25 cm^3) into a clean 100 cm^3 beaker followed by adding one t.t. of distilled water.

Switch on the pH meter and adjust the temperature knob to the lab temperature. Connect the combined electrode (i.e. glass and calomel electrode) to the pH meter. Wash the electrode with distilled water and immerse the electrode into the acid taken in a beaker. Note down the pH of the acid. Add 0.5 cm^3 of sodium hydroxide solution at a time taken in a micro burette in to the acid, stir the solution using magnetic stirrer and measure the pH of the acid and continue the addition of sodium hydroxide until the sudden raise in the pH. Then take few more readings and plot the graph, pH against the volume of sodium hydroxide added (graph No1). Take the difference in the pH of acid (ΔpH) and the volume of sodium hydroxide added (ΔV). Plot the graph, $\Delta\text{pH}/\Delta V$ against the volume of sodium hydroxide added (graph No.2) and from the graph 2, determine the equivalence point of the acid. Determine the half equivalence point that gives the pK_a value of weak acid by the graph 1.

Note: **Switch off the instrument after experiment is completed.**

RESULTS:

1. Equivalence Point = _____
2. Half E.P = _____
3. pH at E.P = _____
4. pH at H.E.P. = _____
5. pK_a of the given acid = _____

OBSERVATIONS AND CALCULATIONS

| Trial No. | I | II | III | Average time in sec |
|---------------------------------------|---|----|-----|---------------------|
| Time taken by water in sec (t_w) | | | | |
| Time taken by liquid in sec (t_L) | | | | |

$$\eta_L = \frac{\eta_w t_L d_L}{t_w d_w}$$

Note: lab temperature=.....°C

Where t_w = time taken by water = _____ sec

t_L = time taken by liquid = _____ sec

d_w = density of water = _____ g/cc

d_L = density of liquid = _____ g/cc

η_w = **co – efficient of viscosity of water** = _____ **milli poise**

η_L = co – efficient of viscosity of liquid = _____ milli poise

$$\eta_L = \frac{\eta_w t_L d_L}{t_w d_w} = \text{_____ milli poise}$$

$\therefore \eta_L$ = **co – efficient of viscosity of liquid** = _____ **milli poise**

EXPERIMENT NO.5Date:

| | | | | | | | |
|--|--|--|--|--|--|--|--|
| | | | | | | | |
|--|--|--|--|--|--|--|--|

**DETERMINATION OF VISCOSITY COEFFICIENT OF GIVEN LIQUID
USING OSTWALD'S VISCOMETER**

THEORY: Viscosity arises due to internal friction between moving layers of molecules. A liquid flowing through a cylindrical tube of uniform diameter is expected to move in the form of molecular layers. A layer close to the surface is almost stationary while that at the axis of tube moves faster than any other intermediate layer. A slow moving layer exerts a drag or friction on its nearest moving layer backwards. This property of a liquid by which it retards or opposes motion of a liquid between the layers is called viscosity. The coefficient of viscosity is defined as the tangential force per unit area required to maintain a unit velocity gradient between any two successive layers of a liquid situated unit distance apart. The coefficient of viscosity of a liquid is given by Poiseuille's formula,

$$\eta = \frac{\pi P r^4 t}{8 V l}, \text{ Where } P = h d g$$

$$\eta = \frac{\pi h d g r^4 t}{8 V l}$$

Where η is the coefficient of viscosity of the liquid, V = volume of the liquid, r = radius of the tube, l = length of the tube, P is the pressure difference between the two ends of the tube, h = height of the capillary tube, d = density of the liquid, and g = gravitational force. If equal volumes of two different liquids are allowed to flow through the same tube under identical

conditions, then,

$$\frac{\eta_w}{\eta_L} = \frac{t_w d_w}{t_L d_L}$$

Where, η_w = co-efficient of viscosity of the standard liquid (water), and η_L = co-efficient of viscosity of the given liquid which can be determined. t_w = time taken by the standard liquid (Water), t_L = time taken by the given liquid, d_w = density of standard liquid (water), d_L = density of given liquid. Viscosity is a temperature dependent property and is inversely proportional to the temperature (T), i.e. $\eta \propto 1/T$. Hence the measurements are carried out at constant temperature.




Fig. Viscometer

PROCEDURE:

1. Clean the given viscometer with distilled water and rinse with acetone. Dry it in an electric oven and then fix the viscometer vertically to a stand.
2. Transfer a known volume (say 15 cm³) of distilled water taken in a burette into the viscometer. Suck the distilled water above the upper mark 'A' in the viscometer and allow it to flow freely through a capillary. Start the stop clock when the water level just crosses the upper mark 'A' and stop it when the water level just crosses the lower mark 'B' in the viscometer. Note down the time of flow of water in seconds (t_w) and repeat the experiment to get agreeing values.
3. Pour out the distilled water and rinse the viscometer with acetone again. Dry it and fix the viscometer vertically to a stand.
4. Transfer the known volume (say 15 cm³) of the given organic liquid into a viscometer. Note down the time of flow of organic liquid (t_L) in seconds as similar to the water and repeat the experiment to get agreeing values.
5. Determine the viscosity co-efficient of the given organic liquid using the relation,

$$\eta_L = \frac{\eta_w t_L d_L}{t_w d_w} \text{ and the given values of } \eta_w, d_w \text{ and } d_L.$$

RESULT: $\therefore \eta_L$ = co – efficient of viscosity of liquid = _____
milli poise

Note: Transfer the organic liquid into the bottle after the experiment is completed and rinse the viscometer with acetone and keep it in an electric oven.

PART-B
VOLUMETRIC
ANALYSIS

OBSERVATIONS AND CALCULATIONS:

1. Preparation of standard solution of disodium salt of EDTA

Weight of weighing bottle + EDTA salt = _____ g.

Weight of empty weighing bottle = _____ g.

∴ Weight of EDTA salt taken = _____ g.

$$\text{Molarity of EDTA} = \frac{\text{Weight of EDTA salt}}{\text{Molecular weight of EDTA}} \times 4$$

$$M_{\text{EDTA}} = \frac{\quad \times 4}{372.24} = \quad \text{(a) M}$$

2. Determination of total hardness of water.

Burette : Standard EDTA solution

Conical flask : 25 cm³ of hard water + 5cm³ of NH₄OH – NH₄Cl buffer .

Indicator : Eriochrome black – T

End point : Wine red to clear blue colour

| Burette readings(cm ³) | I | II | III |
|------------------------------------|---|----|-----|
| Final readings | | | |
| Initial readings | | | |
| Volume of EDTA run down | | | |



∴ Volume of EDTA required V= _____ cm³

EXPERIMENT NO. 1Date:

| | | | | | | | |
|--|--|--|--|--|--|--|--|
| | | | | | | | |
|--|--|--|--|--|--|--|--|

DETERMINATION OF TOTAL HARDNESS OF WATER.**THEORY:**

Hardness of water is due to the presence of dissolved salts of calcium and magnesium. The total hardness of water is determined using ethylene diamine tetra acetic acid (EDTA) which form complex with Ca^{2+} and Mg^{2+} ions. The EDTA molecule has two easily replaceable hydrogen atoms and hence the EDTA form complex with metal ion present in water. Reaction can be carried out quantitatively at a pH of 10 using Eriochrome black-T as indicator. Since the reaction involves the liberation of H^+ ions, a buffer mixture has to be used to maintain a pH of 10. The buffer mixture used in the titration is ammonium hydroxide and ammonium chloride. The total hardness of the water is usually expressed in terms of ppm (ie, parts per million) of CaCO_3 . EDTA is sparingly soluble in water and hence disodium salt of EDTA is used as reagent for preparing the solution. The structures of EDTA (H_4Y) and disodium salt of EDTA ($\text{Na}_2\text{H}_2\text{Y}$) are given below.



NOTE: Molecular weight of $\text{CaCO}_3 = 100$

1000 cm^3 of 1M of EDTA = one g molecular weight of CaCO_3 (100)

$$\therefore V \text{ cm}^3 \text{ of (a) M EDTA} = X = \frac{(a) \times V \times 100}{1000} =$$

$$X = \text{_____g of CaCO}_3$$

$$\therefore \text{Weight of CaCO}_3 \text{ present in } 25 \text{ cm}^3 \text{ of hard water} = X = \text{_____g}$$

\therefore Weight of CaCO_3 present in 25 cm^3

$$\text{of hard water in terms of ppm} = Y = \frac{X \times 10^6}{25}$$

$$Y = \frac{\text{_____} \times 10^6}{25} = \text{_____ppm}$$

Thus total hardness of water in terms of ppm = Y = _____ppm

Note: Parts per million: ppm is the parts of CaCO_3 equivalent hardness present per million **or** 10^6 parts of water.

ie: 500 ppm = 500 part of CaCO_3 equivalent hardness present per million **or** 10^6 parts of water.

PROCEDURE:


I. Preparation of standard solution of disodium salt of EDTA:

Weigh given disodium salt of EDTA accurately using an electronic balance. Transfer the crystals of EDTA into a clean 250 cm³ standard flask using a funnel. To this add **one** test tube of ammonium hydroxide and little distilled water and shake well to dissolve the crystals. Make up the solution to the mark by adding distilled water and mix well for uniform concentration. Calculate the molarity of EDTA solution.

II. Determination of total hardness of water:

Pipette out 25cm³ of the given hard water sample into a clean conical flask, to this add 5 cm³ of ammonia-ammonium chloride buffer solution, and 2 drops of Eriochrome black-T indicator. Titrate this mixture against standard EDTA solution taken in a burette until the wine red color changes to clear blue. Note down the volume of EDTA added and repeat the experiment for concordant values. Calculate the total hardness of water in terms of ppm.

REACTIONS:



RESULT: The total hardness of the given water sample is.....ppm of CaCO₃.

OBSERVATIONS AND CALCULATIONS:

I . Preparation of standard solution of disodium salt of EDTA

Weight of weighing bottle + EDTA salt = _____ g.

Weight of empty weighing bottle = _____ g.

∴ Weight of EDTA salt taken = _____ g.

$$\text{Molarity of EDTA} = \frac{\text{Weight of EDTA salt}}{\text{Molecular weight of EDTA}} \times 4$$

$$M_{\text{EDTA}} = \frac{\quad \times 4}{372.24} = \quad \text{(a) M}$$

II. Estimation of CaO in cement solution

Burette : Standard EDTA solution

Conical flask : 25 cm³ of cement solution + 5cm³ of 1:1 glycerol
+ 5cm³ of di-ethylamine + 10 cm³ of 4N NaOH .

Indicator : Pattern - Reeder's indicator

End point : Wine red to purple blue color

| Burette readings | I | II | III |
|--|---|----|-----|
| Final readings | | | |
| Initial readings | | | |
| Volume of EDTA run down in cm ³ | | | |

∴ Volume of EDTA required = V= _____ cm³

EXPERIMENT NO. 2Date:

| | | | | | | | |
|--|--|--|--|--|--|--|--|
| | | | | | | | |
|--|--|--|--|--|--|--|--|

DETERMINATION OF CALCIUM OXIDE IN THE GIVEN SAMPLE OF CEMENT SOLUTION BY RAPID EDTA METHOD**THEORY:**

Cement contains oxides of calcium, aluminum, magnesium, iron and small amount of silica. Calcium oxide is a prime constituent of cement and the general composition of portland cement is given below.

CaO = 60 – 66%, SiO₂ = 17 – 25%, Al₂O₃ = 3 – 8%, Fe₂O₃ = 2 – 6%, MgO = 0.1- 5.5%, SO₃ = 1 – 3 % and Na₂O & K₂O = 0.5 – 1.5%

In the estimation of calcium in cement, the given cement sample is treated with concentrated hydrochloric acid. The insoluble silica is filtered off and the filtrate which contains calcium ions is titrated against EDTA at a pH of 12 – 14 using Pattern-Reeder's indicator. The latter is used in this titration because Eriochrome black – T indicator forms a very weak complex with calcium ions. Diethylamine is added to maintain a pH of about 12.5 and glycerol is added to get a sharp end point. This method facilitates the determination of calcium only although Mg⁺² ions are present in the cement solution which can readily form complexes with EDTA (i.e., Mg⁺² ions do not react with EDTA in the above conditions)

Procedure:**I. PREPARATION OF STANDARD SOLUTION OF DISODIUM SALT OF EDTA:**

Weigh the given disodium salt of EDTA accurately using an electronic balance. Transfer the crystals of EDTA into a clean 250cm³ standard flask using a funnel. Dissolve the crystals by adding **one** t.t. of ammonium hydroxide and little distilled water. Make up the solution to the mark by adding distilled water and mix well for uniform concentration. Calculate the molarity of EDTA solution.

NOTE: Molecular weight of CaO = 56.08

1000 cm³ of 1M of EDTA = one gm molecular weight of CaO (56.08)

$$\therefore V \text{ cm}^3 \text{ of (a) M EDTA} = X = \frac{(a) \times V \times 56.08}{1000} = \frac{x}{1000} \times 56.08$$

$$X = \text{_____ g of CaO}$$

$$\therefore \text{Weight of CaO present in 25 cm}^3 \text{ of cement solution} = X = \text{_____g}$$

$$\therefore \text{Weight of CaO present in 1000cm}^3 \text{ of cement solution} = Y = \frac{X \times 1000}{25}$$

$$Y = \text{_____}$$

$$Y = \text{_____g}$$



Result:

Amount of CaO in the 1000 cm³ cement solution= _____g

II. Estimation of calcium oxide:

Pipette out 25cm³ of the given cement solution into a clean conical flask. Add 5cm³ of 1:1 glycerol, 5cm³ of diethyl amine and 10cm³ of 4N sodium hydroxide. Add 3-4 drops of Pattern- Reeder's indicator and titrate this against EDTA solution taken in a burette until wine red color changes to purple blue. Note down the volume of EDTA added and repeat the experiment for concordant values. Calculate the percentage of calcium oxide in the given cement solution.

REACTIONS:



RESULT: The amount of CaO present in the given cement solution is

..... g

OBSERVATIONS AND CALCULATIONS:

I. Preparation of brass solution:

Weight of brass taken is =g (W)

II. Estimation of copper:

Burette : Standard hypo(sodium thiosulphate) solution

Conical flask : 25 cm³ brass solution + NH₄OH drop wise till a bluish white ppt. is formed. Dissolve the ppt. by adding dilute acetic acid + $\frac{1}{4}$ th t.t. of acetic acid in excess + One t.t. of 10% KI + shake well + one t.t. of distilled water.

Indicator : 2 cm³ of freshly prepared starch.

End point : Disappearance of blue color

| Burette readings | I | II | III |
|--|---|----|-----|
| Final readings | | | |
| Initial readings | | | |
| Volume of hypo run down in cm ³ . | | | |

∴ Volume of hypo required = V= _____ cm³

NOTE: 1. Normality of hypo = 0.05 N

 2. One gram equivalent weight of copper = 63.54

EXPERIMENT NO. 3**Date:**

| | | | | | | | |
|--|--|--|--|--|--|--|--|
| | | | | | | | |
|--|--|--|--|--|--|--|--|

**DETERMINATION OF PERCENTAGE OF COPPER IN BRASS USING
STANDARD SODIUM THIOSULPHATE SOLUTION.****THEORY :**

The chief constituents of brass alloy are copper and zinc. It also contains small quantities of tin, lead and iron. The percentage composition of a typical brass alloy is given below.

Cu = 50 – 90 %, Zn = 20 – 40 % , Sn = 0 - 6% , Pb = 0 - 2% and Fe = 0 – 1% .

A solution of brass is made by dissolving the sample in minimum amount of nitric acid and excess of nitric acid destroyed by boiling with urea. The solution is neutralized and made in to slight acidic medium. Potassium iodide solution is added and the cupric ion present in the brass solution oxidise potassium iodide to iodine. The liberated iodine is titrated against sodium thiosulphate using starch as indicator. The volume of sodium thiosulphate consumed is measure of the amount of copper present in the solution and sodium thiosulphate which is commonly known as hypo.

PROCEDURE:**I. Preparation of Brass solution:**

Weigh the given brass pieces accurately weigh and transfer it into a clean conical flask. Add one t.t. of 1:1 nitric acid. Warm gently to dissolve the brass piece completely and add two t.t of distilled water. Add pinch of urea and continue the heating until the reddish brown fumes of oxides of nitrogen are expelled completely. Cool the solution to room temperature and transfer the solution of brass in to a clean 250 cm³ standard flask. Make up the solution to the mark by adding distilled water and mix well for uniform concentration.

∴ 1000 cm³ of 1 N of hypo = one g equivalent weight of copper (63.54)

$$\therefore V \text{ cm}^3 \text{ of } 0.05 \text{ N hypo} = X = \frac{V \times 0.05 \text{ N} \times 63.54}{1000} = \text{-----}$$

$$X = \text{-----} \text{ g of copper}$$

∴ Weight of copper present in 25 cm³ of brass solution = X = -----g

$$\therefore \text{Weight of copper present in } 250 \text{ cm}^3 \text{ of brass solution} = Y = \frac{X \times 250}{25}$$

$$Y = \text{-----} \text{ g.}$$

∴ Percentage of copper in the

$$\text{given brass sample} = \frac{Y \times 100}{\text{Weight of brass taken (W)}}$$





$$= \text{-----}$$

∴ Percentage of copper in the given brass sample =-----%

II Estimation of Copper:

Pipette out 25 cm³ of brass solution into a clean conical flask and add dilute ammonium hydroxide drop by drop until bluish white precipitate is formed. Dissolve the ppt. by adding dilute acetic acid drop by drop and then add 1/4th t.t. of acetic acid in excess. Add one t.t. of 10% potassium iodide solution and shake well. The liberated iodine titrated against standard sodium thiosulphate solution taken in a burette using starch as indicator towards the end point. (add starch when brown solution changes to light yellow). Disappearance of blue colour is the end point. Note down the volume of hypo added and repeat the experiment for concordant values. Calculate the percentage of copper in the given sample of brass.

REACTIONS:



RESULT:

The percentage of Cu present in the given brass sample is.....%

OBSERVATIONS AND CALCULATIONS:

I. Preparation of standard potassium dichromate solution :

Weight of weighing bottle + $K_2Cr_2O_7$ crystals = _____ g.

Weight of empty weighing bottle = _____ g.

∴ Weight of $K_2Cr_2O_7$ crystals taken = _____ g.

$$\text{Normality of } K_2Cr_2O_7 = \frac{\text{Weight of } K_2Cr_2O_7 \times 4}{\text{Equivalent weight of } K_2Cr_2O_7}$$

$$N_{K_2Cr_2O_7} = \frac{x \times 4}{49} = \dots\dots\dots N(a)$$

II. Estimation of iron:

Burette : Standard $K_2Cr_2O_7$ solution

Conical flask : 25 cm³ of haematite ore solution + ¼ t.t concentrated HCl + heat + SnCl₂ drop wise till yellow color changes to colorless + 5-6 drops of SnCl₂ in excess + cool the solution to room temperature + 10 cm³ of saturated HgCl₂ + shake well + one t.t. of distilled H₂O

Indicator : potassium ferri cyanide as external indicator

End point : Failure to observe blue color; when a drop of reacting mixture is brought in contact with the indicator placed on paraffine paper

| Burette readings | pilot | I | II | III |
|--|-------|---|----|-----|
| Final readings | | | | |
| Initial readings | | | | |
| Volume of $K_2Cr_2O_7$ run down in cm ³ s | | | | |

∴ Volume of $K_2Cr_2O_7$ required = V = _____ cm³

EXPERIMENT NO. 4**Date:**

| | | | | | | | |
|--|--|--|--|--|--|--|--|
| | | | | | | | |
|--|--|--|--|--|--|--|--|

DETERMINATION OF IRON IN A GIVEN SAMPLE OF HAEMATITE ORE SOLUTION**THEORY:**

Haematite is an important ore of iron containing mainly ferric oxide (Fe_2O_3) and small amount of silica (SiO_2). A known weight of the ore is digested with hydrochloric acid and the insoluble residue mainly containing silica is removed by filtration. The filtrate is diluted to 250 cm^3 and is supplied for the estimation. Iron is present in the solution as Fe^{3+} (ferric) ions and it is reduced to Fe^{2+} (ferrous) ions using stannous chloride (SnCl_2) in hot condition in presence of concentrated HCl . The excess of stannous chloride is oxidized by adding mercuric chloride solution. The ferrous ions in the resulting solution is titrated against standard potassium dichromate using potassium ferricyanide as external indicator. Potassium ferricyanide produces an intense deep blue color with ferrous ions due to the formation of ferriferrocyanide complex.

PROCEDURE:**a) Preparation of standard solution of potassium dichromate:**

Weigh given potassium dichromate crystals accurately using an electronic balance. Transfer the crystals of potassium dichromate into a clean 250 cm^3 standard flask using a funnel. Dissolve the crystals by adding distilled water and make up the solution to the mark by adding distilled water. Mix well for uniform concentration and calculate the normality of potassium dichromate solution.

NOTE: **Equivalent weight of iron** = **55.85 g**

1000 cm³ of 1 N of K₂Cr₂O₇ = one gm equivalent weight of iron (55.85g)

$$\therefore \text{V cm}^3 \text{ of (a) N}_{\text{K}_2\text{Cr}_2\text{O}_7} = \text{X} = \frac{(\text{a}) \times \text{V} \times 55.85}{1000}$$

$$X = \frac{\quad}{1000} = \quad gm \text{ of iron}$$

∴ Weight of iron present in 25 cm³ of haematite solution = X = _____ g.




\therefore Weight of iron present in 1000 cm^3 of

$$\text{Haematite ore solution} = Y = \frac{X \times 1000}{25} =$$

Y = _____g.

b) Estimation of iron:

Pipette out 25cm³ of the haematite ore solution into a clean conical flask. Add ¼ t.t of concentrated hydrochloric acid and heat the solution nearly to boiling. Add stannous chloride drop by drop taken in a separate burette till yellow color solution becomes colorless and then add 6-8 drops of stannous chloride in excess. Cool the solution to lab temperature and add 10cm³ of saturated mercuric chloride at once, shake well and a silky white precipitate is formed. (if no ppt or black ppt is formed, reject and repeat the experiment).add one t.t. of distilled water and titrate this against standard potassium dichromate solution taken in a burette, using potassium ferricyanide as external indicator. the end point is failure to observe blue colour when one drop of the titrated mixture is brought in contact with a drop of the indicator placed on paraffin paper. note down the volume of potassium dichromate added and repeat the experiment for concordant values. calculate the amount of iron in the given sample of haematite ore solution.

REACTIONS:**RESULT:**

The amount of iron in the given Haematite Ore solution is..... g

OBSERVATIONS AND CALCULATIONS:

I. Preparation of Mohr's salt or ferrous ammonium sulphate (FAS) solution:

Weight of weighing bottle + FAS crystals = _____ g.

Weight of empty weighing bottle = _____ g.

∴ Weight of FAS crystals taken = _____ g.

$$\text{Normality of FAS} = \frac{\text{Weight of FAS} \times 4}{\text{Equivalent weight of FAS}}$$

$$N_{\text{FAS}} = \frac{\text{_____} \times 4}{392} = \text{_____ (a) N}$$

II. Estimation of Chemical Oxygen Demand(back titration)

Burette : Standard FAS solution

Conical flask : 25 cm³ of waste water +25 cm³ of standard K₂Cr₂O₇ solution + one tt of 1:1 H₂SO₄ .

Indicator : 4 – 6 drops of ferroin indicator

End point : Blue green to reddish brown color

| | | | | |
|--|---|----|-----|--|
| Burette readings | I | II | III | |
| Final readings | | | | |
| Initial readings | | | | |
| Volume of FAS run down in cm ³ | | | | |

∴ Volume of FAS solution required = V₁ = _____ cm³

EXPERIMENT NO. 5Date:

| | | | | | | | |
|--|--|--|--|--|--|--|--|
| | | | | | | | |
|--|--|--|--|--|--|--|--|

**DETERMINATION OF CHEMICAL OXYGEN DEMAND (COD) OF THE
GIVEN INDUSTRIAL WASTE WATER SAMPLE****THEORY:**

The chemical oxygen demand (COD) test is extensively employed for the measurement of the pollution strength of industrial wastes. Chemical oxygen demand is a measure of the total quantity of oxygen required for oxidation of organic compounds of wastes to CO_2 and water by a strong oxidizing agent. This parameter is particularly valuable in surveys designed to determine the control losses to sewer systems. Results may be obtained within a relatively short time and measures taken to correct error on the day they occur.

Waste water contains organic impurities which include straight chain aliphatic and aromatic compounds, such as alcohols, acids, amines, pyridine and other oxidisable materials. They oxidize more effectively when silver sulphate is added as a catalyst. But silver sulphate reacts with chlorides in the waste water to form precipitates which are oxidized partially by this procedure. This difficulty is overcome by adding mercuric sulphate to the sample.

PROCEDURE:**I. PREPARATION OF STANDARD MOHR'S SALT OR FERROUS
AMMONIUM SULPHATE (FAS) SOLUTION:**

Weigh given Mohr's salt crystals accurately using an electronic balance and transfer the crystals into a clean 250 cm^3 standard flask. Dissolve the crystals by adding two t.t. of dilute sulphuric acid and little distilled water. and make it up to the mark by adding distilled water. Mix well for uniform concentration and calculate the normality of Mohr's salt solution.

IV Blank titration or (FAS > < K₂Cr₂O₇)

Burette : Standard FAS solution
Conical flask : 25 cm³ of K₂Cr₂O₇ solution + 1tt of 1:1 H₂SO₄
Indicator : 4 – 6 drops of Ferrion indicator
End Point : Blue green to reddish brown color

∴ Volume of FAS solution required for blank titration
= **V₂** = Final B.R - Initial B.R=.....

$$V_2 = \text{_____ cm}^3$$

∴ Volume of K₂Cr₂O₇ required for the oxidation
of organic compounds of waste water = **V₂ - V₁** = _____ = _____ cm³

1000 cm³ of 1 N of FAS = 1000 cm³ of 1N K₂Cr₂O₇ solution = 8 g of oxygen

$$V_2 - V_1 \text{ cm}^3 \text{ of (a) N FAS} = X = \frac{8 \times (\text{a}) \times V_2 - V_1}{1000} =$$

$$X = \frac{8 \times \text{_____}}{1000}$$

$$X = \text{_____ g of oxygen}$$

∴ Amount of O₂ required to oxidise organic
compounds present in 25 cm³ of waste water = X = _____ g.

∴ Amount of O₂ required to oxidize organic
compounds present in 1000cm³ of waste water = Y = $\frac{X \times 1000}{25}$ = _____

$$Y = \text{_____ g.}$$

∴ COD of waste water = Y x 1000 = _____ x 1000 mg

∴ COD of waste water = _____ mg of oxygen


II. Determination of COD of waste water:

Pipette out 25 cm³ of the industrial waste water and add exactly 25cm³ of standard potassium dichromate solution (0.1N) into a clean conical flask. Add one tt of 1:1 sulphuric acid and 2 to 4 drop of ferroin indicator and titrate against standard Mohr's salt solution taken in a burette until blue-green color changes to reddish-brown. Note down the volume of Mohr's salt solution added for unreacted potassium dichromate solution (V₁) and repeat the experiment to get concordant values. Using blank titre value determines the volume of potassium dichromate required to oxidize organic compounds and calculate the COD of the given industrial waste water.

III. Blank titration

Add exactly 25cm³ of standard potassium dichromate solution into a clean conical flask and one tt of 1:1 sulphuric acid . Add 2–4 drops of ferrion indicator and titrate this against Mohr's salt solution taken in a burette until blue–green colour changes to reddish brown. Note down the volume of Mohr's salt solution required for blank titration (V₂).

REACTIONS:



RESULT:

The COD of the given waste water sample containsmg of oxygen

Determination of chloride content of water (Iodometric method)

a. Principle: Chlorine will liberate free iodine from potassium iodide (KI) solutions at pH 8 or less. The liberated iodine is titrated with a standard solution of sodium thiosulfate ($\text{Na}_2\text{S}_2\text{O}_3$) with starch as the indicator. Titrate at pH 3 to 4 because the reaction is not stoichiometric at neutral pH due to partial oxidation of thiosulfate to sulfate.

b. Interference: Oxidized forms of manganese and other oxidizing agents interfere. Reducing agents such as organic sulfides also interfere. Although the neutral titration minimizes the interfering effect of ferric and nitrite ions, the acid titration is preferred because some forms of combined chlorine do not react at pH 7. Use only acetic acid for the acid titration; sulfuric acid (H_2SO_4) will increase interferences; *never use hydrochloric acid (HCl)*. See Section A.3 for discussion of other interferences.

2. Apparatus

For a description of the amperometric end-point detection apparatus and a discussion of its use, see D.2a.

3. Reagents

- a. Standard phenylarsine oxide solution, 0.005 64N:* Dissolve approximately 0.8 g phenylarsine oxide powder in 150 mL 0.3N NaOH solution. After settling, decant 110 mL into 800 mL distilled water and mix thoroughly. Bring to pH 6 to 7 with 6N HCl and dilute to 950 mL with distilled water.

CAUTION: *Severe poison, cancer suspect agent.*

Standardization—Accurately measure 5 to 10 mL freshly standardized 0.0282N iodine solution into a flask and add 1 mL KI solution. Titrate with phenylarsine oxide solution, using the amperometric end point (Method D) or starch solution (see B.2e) as an indicator. Adjust to 0.005 64N and recheck against the standard iodine solution; 1.00 mL = 200 μg available chlorine. (CAUTION: *Toxic—take care to avoid ingestion.*)

Standard sodium thiosulfate solution.

- b. *Standard sodium thiosulfate solution, 0.005 64N:* Prepare by diluting 0.1N $\text{Na}_2\text{S}_2\text{O}_3$. For maximum stability of the dilute solution, prepare by diluting an aged 0.1N solution with freshly boiled distilled water (to minimize bacterial action) and add 4 g $\text{Na}_4\text{B}_4\text{O}_7/\text{L}$. To inhibit mold formation optionally add either 10 mg HgI_2 or 2 drops toluene per liter of solution. Standardize daily as directed in B.2c using 0.005 64N $\text{K}_2\text{Cr}_2\text{O}_7$ or iodate solution. Use sufficient volume of sample so that the final dilution does not exceed 1 + 2. Use an automatic buret of a type in which rubber does not come in contact with the solution. 1.00 mL = 200 μg available chlorine.
- c. *Potassium iodide, KI, crystals.*
- d. *Acetate buffer solution, pH 4.0:* Dissolve 146 g anhydrous $\text{NaC}_2\text{H}_3\text{O}_2$, or 243 g $\text{NaC}_2\text{H}_3\text{O}_2 \cdot 3\text{H}_2\text{O}$, in 400 mL distilled water, add 480 g conc acetic acid, and dilute to 1 L with chlorine-demand-free water.
- e. *Standard arsenite solution, 0.1N:* Accurately weigh a stoppered weighing bottle containing approximately 4.95 g arsenic trioxide, As_2O_3 . Transfer without loss to a 1-L volumetric flask and again weigh bottle. Do not attempt to brush out adhering oxide. Moisten As_2O_3 with water and add 15 g NaOH and 100 mL distilled water. Swirl flask contents gently to dissolve. Dilute to 250 mL with distilled water and saturate with CO_2 , thus converting all NaOH to NaHCO_3 . Dilute to mark, stopper, and mix thoroughly. This solution will preserve its titer almost indefinitely. (CAUTION: *Severe poison. Cancer suspect agent.*)

$$\text{Normality} = \frac{\text{g As}_2\text{O}_3}{49.455}$$

Standard iodine solution, 0.1N: Dissolve 40 g KI in 25 mL chlorine-demand-free water, add 13 g resublimed iodine, and stir until dissolved. Transfer to a 1-L volumetric flask and dilute to mark.

Standardization—Accurately measure 40 to 50 mL 0.1N arsenite solution into a flask and titrate with 0.1N iodine solution, using starch solution as indicator. To obtain accurate results, insure that the solution is saturated with CO_2 at end of titration by passing current of CO_2 through solution for a few minutes just before end point is reached, or add a few drops of HCl to liberate sufficient CO_2 to saturate solution. Alternatively standardize against $\text{Na}_2\text{S}_2\text{O}_3$;

Optionally, prepare 0.1000*N* iodine solution directly as a standard solution by weighing 12.69 g primary standard resublimed iodine. Because I₂ may be volatilized and lose from both solid and solution, transfer the solid immediately to KI as Specified above

Never let solution stand in open containers for extended periods.

*Standard iodine titrant, 0.0282*N**: Dissolve 25 g KI in a little distilled water in a 1-L volumetric flask, add correct amount of 0.1*N* iodine solution exactly standardized to yield a 0.0282*N* solution, and dilute to 1 L with chlorine-demand-free water. For accurate work, standardize daily according to directions in ¶ 3g above, using 5 to 10 mL of arsenite or Na₂S₂O₃ solution. Store in amber bottles or in the dark; protect solution from direct sunlight at all times and keep from all contact with rubber.

- a. *Starch indicator*: See B.2e.
- b. *Standard iodate titrant, 0.005 64*N**: Dissolve 201.2 mg primary standard grade KIO₃, dried for 1 h at 103°C, or 183.3 mg primary standard anhydrous potassium bi-iodate in distilled water and dilute to 1 L.
- c. *Phosphoric acid solution, H₃PO₄, 1 + 9*.
- d. *Phosphoric acid-sulfamic acid solution*: Dissolve 20 g NH₂SO₃H in 1 L 1 + 9 phosphoric acid.
- e. *Chlorine-demand-free water*: Prepare chlorine-demand-free water from good-quality distilled or deionized water by adding sufficient chlorine to give 5 mg/L free chlorine. After standing 2 d this solution should contain at least 2 mg/L free chlorine; if not, discard and obtain better-quality water. Remove remaining free chlorine by placing container in sunlight or irradiating with an ultraviolet lamp. After several hours take sample, add KI, and measure total chlorine with a colorimetric method using a nessler tube to increase sensitivity. Do not use before last trace of free and combined chlorine has been removed.

Distilled water commonly contains ammonia and also may contain reducing agents. Collect good-quality distilled or deionized water in a sealed container from which water can be drawn by gravity. To the air inlet of the container add an H₂SO₄ trap consisting of a large test tube half filled with 1 + 1 H₂SO₄ connected in series with a similar but empty test tube. Fit both test tubes with stoppers and inlet tubes terminating near the bottom of the tubes and outlet tubes terminating near the top of the tubes. Connect outlet tube of trap containing H₂SO₄ to the distilled water container, connect inlet tube to outlet of empty test tube. The empty test tube will prevent discharge to the atmosphere of

H₂SO₄ due to temperature-induced pressure changes. Stored in such a container, chlorine-demand-free water is stable for several weeks unless bacterial growth occurs.

. Preparation for titration:

- 1) Volume of sample—For chlorine concentration of 10 mg/L or less, titrate 200 mL. For greater chlorine concentrations, use proportionately less sample and dilute to 200 mL with chlorine-demand-free water. Use a sample of such size that not more than 10 mL phenylarsine oxide solution is required.
- 2) Preparation for titration—Measure 5 mL 0.005 64*N* phenylarsine oxide or thiosulfate for chlorine concentrations from 2 to 5 mg/L, and 10 mL for concentrations of 5 to 10 mg/L, into a flask or casserole for titration with standard iodine or iodate. Start stirring. For titration by amperometry or standard iodine,

also add excess KI (approximately 1 g) and 4 mL acetate buffer solution or enough to reduce the pH to between 3.5 and 4.2.

b. Titration: Use one of the following:

- 1) Amperometric titration—Add 0.0282*N* iodine titrant in small increments from a 1-mL buret or pipet. Observe meter needle response as iodine is added: the pointer remains practically stationary until the end point is approached, whereupon each iodine increment causes a temporary deflection of the microammeter, with the pointer dropping back to its original position. Stop titration at end point when a small increment of iodine titrant gives a definite pointer deflection upscale and the pointer does not return promptly to its original position. Record volume of iodine titrant used to reach end point.
- 2) Colorimetric (iodine) titration—Add 1 mL starch solution and titrate with 0.0282*N* iodine to the first appearance of blue color that persists after complete mixing.
- 3) Colorimetric (iodate) titration—To suitable flask or casserole add 200 mL chlorine-demand-free water and add, with agitation, the required volume of reductant, an excess of KI (approximately 0.5 g), 2 mL 10% H₃PO₄ solution, and 1 mL starch solution in the order given, and titrate immediately* with 0.005 64*N* iodate solution to the first appearance of a blue color that persists after complete mixing. Designate volume of iodate solution used as *A*. Repeat procedure, substituting 200 mL sample for the 200 mL chlorine-demand-free water. If sample is colored or turbid, titrate to the first change in color, using for comparison another portion of sample with H₃PO₄ added. Designate this volume of iodate solution as *B*.

5. Calculation

a. Titration with standard iodine:

$$\text{mg Cl as Cl}_2/\text{L} = \frac{(A - 5B) \times 200}{C}$$

where:

A = mL 0.005 64*N* reductant,

B = mL 0.0282 *N* I₂, and

C = mL sample.



|| Jai Sri Gurudev ||

BGS INSTITUTE OF TECHNOLOGY

Department of Engineering Chemistry

LABORATORY RUBRICS

| Programme | Course Code | Subject | Credits | L-T-P-TL | Assessment | | Exam Duration |
|-----------|-------------|---------------------------|---------|----------|------------|-----|---------------|
| | | | | | SEE | CIA | |
| B.E | 18CHEL16/26 | Engineering Chemistry Lab | 02 | 0-0-3-3 | 60 | 40 | 3Hrs |

Maximum Marks: 40

| Continuous Internal Evaluation | Excellent (80%-100%) | Good (80%-60%) | Average (40%-50%) |
|--|--|---|---|
| a. Observation write up and punctuality (05) | Students should write the experiments in the Observation book neatly and attend the labs regularly | Students should write the experiments in the Observation book and attend the labs. | Improper maintenance of observation books and being irregular to the labs. |
| b. Conduction of experiment and utput (10) | Students should conduct the experiments following the given procedure, plot the graph, perform calculation and show the accurate results with S.I unit. | Students should conduct the experiments following the given procedure, plot the graph and perform calculation with average results. | Improper conduction of experiments, graph plotting and results without S.I. unit. |
| c. Viva voce (05) | They should answer all the questions. | If they answer some of the questions. | If they doesn't answer the questions. |
| d. Record write up (10) | They should write records neatly, legibly and with suitable Tabular columns. | They should write records with suitable tabular columns | Improper/poor maintenance of record. |
| e. Internal Test (10) | Students should write the given experiments containing Formula, Tabular column, Nature of the graph, conduct the experiment and show the results with S.I. unit. | Students must write the given experiments, conduct the experiment and show the results. | If the student write the experiment but fails to conduct it. |